

Human RETN (Resistin) Antibody Pair (STJAP00056)

STJAP00056

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This RETN antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56729
Gene Symbol	RETN
Uniprot ID	RETN_HUMAN
Immunogen	Capture immunogen: Recombinant Human RETN protein. Detection Immunogen: Recombinant Human RETN protein.
Immunogen Region	
Specificity	
Immunogen Sequence	