

Human VF (Nicotinamide Phosphoribosyltransferase) Antibody Pair (STJAP00063)

STJAP00063

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This VF antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10135
Gene Symbol	NAMPT
Uniprot ID	NAMPT_HUMAN
Immunogen	Capture immunogen: Recombinant Human VF protein. Detection Immunogen: Recombinant Human VF protein.
Immunogen Region	
Specificity	
Immunogen Sequence	