

Mouse EGF (Pro-Epidermal Growth Factor) Antibody Pair (STJAP00068)

STJAP00068

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Antibody Pairs |
| Short Description | This EGF antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA. |
| Applications | ELISA |
| Host/Source | Rat/Goat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal/Polyclonal |
| Clone ID | |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated/Biotin |
| Purification | Capture Purification: Protein A or G. Detection Purification: Antigen Affinity. |
| Dilution Range | |
| Formulation | Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. |
| Isotype | |
| Storage Instruction | Store at -20°C for long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 13645 |
| Gene Symbol | Egf |
| Uniprot ID | EGF_MOUSE |
| Immunogen | Capture immunogen: Recombinant Mouse EGF protein. Detection Immunogen: Recombinant Mouse EGF protein. |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |