

Mouse MMP-9 (Matrix Metalloproteinase-9) Antibody Pair (STJAP00082)

STJAP00082

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This MMP-9 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rat/Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	17395
Gene Symbol	Mmp9
Uniprot ID	MMP9_MOUSE
Immunogen	Capture immunogen: Recombinant Mouse MMP-9 protein. Detection Immunogen: Recombinant Mouse MMP-9 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	