

Mouse TGF-Beta 1 (Transforming Growth Factor Beta-1 Proprotein Cleaved Into-Latency-Associated Peptide) Antibody Pair (STJAP00083)

STJAP00083

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This TGF-Beta 1 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Chicken
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	21803
Gene Symbol	Tgfb1
Uniprot ID	TGFB1_MOUSE
Immunogen	Capture immunogen: Recombinant Mouse TGF-Beta 1 protein. Detection Immunogen: Recombinant Mouse TGF-Beta 1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	