

Rat CRP (C-Reactive Protein) Antibody Pair (STJAP00088)

STJAP00088

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This CRP antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25419
Gene Symbol	Crp
Uniprot ID	CRP_RAT
Immunogen	Capture immunogen: Recombinant Rat CRP protein. Detection Immunogen: Recombinant Rat CRP protein.
Immunogen Region	
Specificity	
Immunogen Sequence	