

Rat KIM-1 (Hepatitis A Virus Cellular Receptor 1 Homolog) Antibody Pair (STJAP00096)

STJAP00096

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This KIM-1 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Goat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	286934
Gene Symbol	Havcr1
Uniprot ID	HAVR1_RAT
Immunogen	Capture immunogen: Recombinant Rat KIM-1 protein. Detection Immunogen: Recombinant Rat KIM-1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	