

Human TSH (Thyrotropin Subunit Beta) Antibody Pair (STJAP00129)

STJAP00129

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This TSH antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A. Detection Purification: Protein A.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7252
Gene Symbol	TSHB
Uniprot ID	TSHB_HUMAN
Immunogen	Capture immunogen: Native Protein. Detection Immunogen: Native Protein.
Immunogen Region	
Specificity	
Immunogen Sequence	