

Human cTnT/TNNT2 (Troponin T-Cardiac Muscle) Antibody Pair (STJAP00137)

STJAP00137

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This cTnT/TNNT2 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7139
Gene Symbol	TNNT2
Uniprot ID	TNNT2_HUMAN
Immunogen	Capture immunogen: Recombinant Human cTnT/TNNT2 protein. Detection Immunogen: Recombinant Human cTnT/TNNT2 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	