

Human ANXA1 (Annexin A1) Antibody Pair (STJAP00178)

STJAP00178

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This ANXA1 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	301
Gene Symbol	ANXA1
Uniprot ID	ANXA1_HUMAN
Immunogen	Capture immunogen: Recombinant Human ANXA1 protein. Detection Immunogen: Recombinant Human ANXA1 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	