

## Human CXCL1 (Growth-Regulated Alpha Protein) Antibody Pair (STJAP00194)

STJAP00194

### GENERAL INFORMATION

<b>Product Type</b>	Antibody Pairs
<b>Short Description</b>	This CXCL1 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Goat/Goat
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal/Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated/Biotin
<b>Purification</b>	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
<b>Dilution Range</b>	
<b>Formulation</b>	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">2919</a>
<b>Gene Symbol</b>	<a href="#">CXCL1</a>
<b>Uniprot ID</b>	<a href="#">GROA_HUMAN</a>
<b>Immunogen</b>	Capture immunogen: Recombinant Human CXCL1 protein. Detection Immunogen: Recombinant Human CXCL1 protein.
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	