

## Human Flt3L (Fms-Related Tyrosine Kinase 3 Ligand) Antibody Pair (STJAP00197)

STJAP00197

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Antibody Pairs  |
| <b>Short Description</b> | This Flt3L antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA. |
| <b>Applications</b>      | ELISA   |
| <b>Host/Source</b>       | Mouse/Goat  |
| <b>Reactivity</b>        | Human   |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal/Polyclonal  |
| <b>Clone ID</b>            |  |
| <b>Concentration</b>       | 0.5 mg/mL  |
| <b>Conjugation</b>         | Unconjugated/Biotin  |
| <b>Purification</b>        | Capture Purification: Protein A or G.<br>Detection Purification: Antigen Affinity.   |
| <b>Dilution Range</b>      |  |
| <b>Formulation</b>         | Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.<br>Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. |
| <b>Isotype</b>             |  |
| <b>Storage Instruction</b> | Store at -20°C for long term, and avoid freeze-thaw cycles.  |

### TARGET INFORMATION

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            | <a href="#">2323</a>   |
| <b>Gene Symbol</b>        | <a href="#">FLT3LG</a>   |
| <b>Uniprot ID</b>         | <a href="#">FLT3L_HUMAN</a>  |
| <b>Immunogen</b>          | Capture immunogen: Recombinant Human Flt3L protein.<br>Detection Immunogen: Recombinant Human Flt3L protein. |
| <b>Immunogen Region</b>   |  |
| <b>Specificity</b>        |  |
| <b>Immunogen Sequence</b> |  |