

## Human Rac1 (Ras-Related C3 Botulinum Toxin Substrate 1) Antibody Pair (STJAP00215)

STJAP00215

### GENERAL INFORMATION

<b>Product Type</b>	Antibody Pairs
<b>Short Description</b>	This Rac1 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Rabbit/Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal/Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated/Biotin
<b>Purification</b>	Capture Purification: Protein A & Antigen Affinity. Detection Purification: Protein A & Antigen Affinity.
<b>Dilution Range</b>	
<b>Formulation</b>	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">5879</a>
<b>Gene Symbol</b>	<a href="#">RAC1</a>
<b>Uniprot ID</b>	<a href="#">RAC1_HUMAN</a>
<b>Immunogen</b>	Capture immunogen: Recombinant Human Rac1 protein. Detection Immunogen: Recombinant Human Rac1 protein.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	