

Human SCF (Kit Ligand) Antibody Pair (STJAP00218)

STJAP00218

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This SCF antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rabbit/Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A. Detection Purification: Protein A.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4254
Gene Symbol	KITLG
Uniprot ID	SCF_HUMAN
Immunogen	Capture immunogen: Recombinant Human SCF protein. Detection Immunogen: Recombinant Human SCF protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	