

Human RET (Proto-Oncogene Tyrosine-Protein Kinase Receptor Ret) Antibody Pair (STJAP00221)

STJAP00221

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This RET antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5979
Gene Symbol	RET
Uniprot ID	RET_HUMAN
Immunogen	Capture immunogen: Recombinant Human RET protein. Detection Immunogen: Recombinant Human RET protein.
Immunogen Region	
Specificity	
Immunogen Sequence	