

Human FGF21 (Fibroblast Growth Factor 21) Antibody Pair (STJAP00223)

STJAP00223

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This FGF21 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	26291
Gene Symbol	FGF21
Uniprot ID	FGF21_HUMAN
Immunogen	Capture immunogen: Recombinant Human FGF21 protein. Detection Immunogen: Recombinant Human FGF21 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	