

Human sRANKL (sRANKL) Antibody Pair (STJAP00226)

STJAP00226

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This sRANKL antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A. Detection Purification: Protein A.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Capture immunogen: Recombinant Human sRANKL protein. Detection Immunogen: Recombinant Human sRANKL protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	