

Human AXL (Tyrosine-Protein Kinase Receptor Ufo) Antibody Pair (STJAP00231)

STJAP00231

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This AXL antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rabbit/Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A & Antigen Affinity. Detection Purification: Protein A & Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	558
Gene Symbol	AXL
Uniprot ID	UFO_HUMAN
Immunogen	Capture immunogen: Recombinant Human AXL protein. Detection Immunogen: Recombinant Human AXL protein.
Immunogen Region	
Specificity	
Immunogen Sequence	