

Human sTREM-1 (Triggering Receptor Expressed On Myeloid Cells 1) Antibody Pair (STJAP00235)

STJAP00235

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This sTREM-1 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54210
Gene Symbol	TREM1
Uniprot ID	TREM1_HUMAN
Immunogen	Capture immunogen: Recombinant Human sTREM-1 protein. Detection Immunogen: Recombinant Human sTREM-1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	