

Human PDGF-AA (Platelet-Derived Growth Factor Subunit A) Antibody Pair (STJAP00250)

STJAP00250

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This PDGF-AA antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5154
Gene Symbol	PDGFA
Uniprot ID	PDGFA_HUMAN
Immunogen	Capture immunogen: Recombinant Human PDGF-AA protein. Detection Immunogen: Recombinant Human PDGF-AA protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	