

Human CCL18 (C-C Motif Chemokine 18) Antibody Pair (STJAP00255)

STJAP00255

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Antibody Pairs |
| Short Description | This CCL18 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA. |
| Applications | ELISA |
| Host/Source | Rabbit/Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | Polyclonal/Polyclonal |
| Clone ID | |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated/Biotin |
| Purification | Capture Purification: Protein A. Detection Purification: Protein A. |
| Dilution Range | |
| Formulation | Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. |
| Isotype | |
| Storage | Store at -20°C for long term, and avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|--|
| Gene ID | 6362 |
| Gene Symbol | CCL18 |
| Uniprot ID | CCL18_HUMAN |
| Immunogen | Capture immunogen: Recombinant Human CCL18 protein. Detection Immunogen: Recombinant Human CCL18 protein. |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |