

Human CYFRA21-1 (Keratin-Type I Cytoskeletal 19) Antibody Pair (STJAP00262)

STJAP00262

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This CYFRA21-1 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3880
Gene Symbol	KRT19
Uniprot ID	K1C19_HUMAN
Immunogen	Capture immunogen: Recombinant Human CYFRA21-1 protein. Detection Immunogen: Recombinant Human CYFRA21-1 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	