

Mouse ANGPTL3 (Angiopoietin-Related Protein 3) Antibody Pair (STJAP00268)

STJAP00268

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This ANGPTL3 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rat/Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	30924
Gene Symbol	Angptl3
Uniprot ID	ANGL3_MOUSE
Immunogen	Capture immunogen: Recombinant Mouse ANGPTL3 protein. Detection Immunogen: Recombinant Mouse ANGPTL3 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	