

Mouse RETN (Resistin) Antibody Pair (STJAP00271)

STJAP00271

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This RETN antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rat/Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	57264
Gene Symbol	Retn
Uniprot ID	RETN_MOUSE
Immunogen	Capture immunogen: Recombinant Mouse RETN protein. Detection Immunogen: Recombinant Mouse RETN protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	