

Mouse INS (INS) Antibody Pair (STJAP00276)

STJAP00276

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This INS antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	16333 16334
Gene Symbol	Ins1 Ins2
Uniprot ID	INS1_MOUSE INS2_MOUSE
Immunogen	Capture immunogen: Recombinant Mouse INS protein. Detection Immunogen: Recombinant Mouse INS protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	