

Mouse IGFBP-5 (Insulin-Like Growth Factor-Binding Protein 5) Antibody Pair (STJAP00277)

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This IGFBP-5 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	16011
Gene Symbol	Igfbp5
Uniprot ID	IBP5_MOUSE
Immunogen	Capture immunogen: Recombinant Mouse IGFBP-5 protein. Detection Immunogen: Recombinant Mouse IGFBP-5 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	