

Mouse MPIF2 (C-C Motif Chemokine 24) Antibody Pair (STJAP00279)

STJAP00279

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This MPIF2 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rat/Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	56221
Gene Symbol	Ccl24
Uniprot ID	CCL24_MOUSE
Immunogen	Capture immunogen: Recombinant Mouse MPIF2 protein. Detection Immunogen: Recombinant Mouse MPIF2 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	