

Mouse EGFR (Epidermal Growth Factor Receptor) Antibody Pair (STJAP00290)

STJAP00290

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This EGFR antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rabbit/Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A & Antigen Affinity. Detection Purification: Protein A & Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13649
Gene Symbol	Egfr
Uniprot ID	EGFR_MOUSE
Immunogen	Capture immunogen: Recombinant Mouse EGFR protein. Detection Immunogen: Recombinant Mouse EGFR protein.
Immunogen Region	
Specificity	
Immunogen Sequence	