

## Mouse SHH (Sonic Hedgehog Protein) Antibody Pair (STJAP00308)

STJAP00308

### GENERAL INFORMATION

<b>Product Type</b>	Antibody Pairs
<b>Short Description</b>	This SHH antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Rat/Rat
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal/Monoclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated/Biotin
<b>Purification</b>	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
<b>Dilution Range</b>	
<b>Formulation</b>	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	20423
<b>Gene Symbol</b>	Shh
<b>Uniprot ID</b>	SHH_MOUSE
<b>Immunogen</b>	Capture immunogen: Recombinant Mouse SHH protein. Detection Immunogen: Recombinant Mouse SHH protein.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	