

Rat SELL (L-Selectin) Antibody Pair (STJAP00333)

STJAP00333

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This SELL antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Goat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	Sell
Uniprot ID	LYAM1_RAT
Immunogen	Capture immunogen: Recombinant Rat SELL protein. Detection Immunogen: Recombinant Rat SELL protein.
Immunogen Region	
Specificity	
Immunogen Sequence	