

Rat EGF (Pro-Epidermal Growth Factor) Antibody Pair (STJAP00337)

STJAP00337

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This EGF antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Sheep/Sheep
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25313
Gene Symbol	Egf
Uniprot ID	EGF_RAT
Immunogen	Capture immunogen: Recombinant Rat EGF protein. Detection Immunogen: Recombinant Rat EGF protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	