

Rat NGF (Beta-Nerve Growth Factor) Antibody Pair (STJAP00341)

STJAP00341

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This NGF antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Goat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	310738
Gene Symbol	Ngf
Uniprot ID	NGF_RAT
Immunogen	Capture immunogen: Recombinant Rat NGF protein. Detection Immunogen: Recombinant Rat NGF protein.
Immunogen Region	
Specificity	
Immunogen Sequence	