

Rat sCD86 (sCD86) Antibody Pair (STJAP00342)

STJAP00342

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This sCD86 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A. Detection Purification: Protein A.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Capture immunogen: Recombinant Rat sCD86 protein. Detection Immunogen: Recombinant Rat sCD86 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	