

Rat GRO Beta/CXCL2 (C-X-C Motif Chemokine 2) Antibody Pair (STJAP00345)

STJAP00345

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This GRO Beta/CXCL2 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Goat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	114105
Gene Symbol	Cxcl2
Uniprot ID	CXCL2_RAT
Immunogen	Capture immunogen: Recombinant Rat GRO Beta/CXCL2 protein. Detection Immunogen: Recombinant Rat GRO Beta/CXCL2 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	