

Rat IL-12 (IL-12) Antibody Pair (STJAP00349)

STJAP00349

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This IL-12 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rabbit/Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Capture immunogen: Recombinant Rat IL-12 protein. Detection Immunogen: Recombinant Rat IL-12 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	