

Rat SAP (Serum Amyloid P-Component) Antibody Pair (STJAP00359)

STJAP00359

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This SAP antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Goat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	29339
Gene Symbol	Apcs
Uniprot ID	SAMP_RAT
Immunogen	Capture immunogen: Recombinant Rat SAP protein. Detection Immunogen: Recombinant Rat SAP protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	