

Rat HPX (Hemopexin) Antibody Pair (STJAP00363)

STJAP00363

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This HPX antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rabbit/Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A. Detection Purification: Protein A & Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	58917
Gene Symbol	Hpx
Uniprot ID	HEMO_RAT
Immunogen	Capture immunogen: Recombinant Rat HPX protein. Detection Immunogen: Recombinant Rat HPX protein.
Immunogen Region	
Specificity	
Immunogen Sequence	