

Rat VWF (Von Willebrand Factor) Antibody Pair (STJAP00365)

STJAP00365

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This VWF antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Sheep/Sheep
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	Vwf
Uniprot ID	VWF_RAT
Immunogen	Capture immunogen: Native Protein. Detection Immunogen: Native Protein.
Immunogen Region	
Specificity	
Immunogen Sequence	