

## Human ANGPTL3 (Angiopoietin-Related Protein 3) Antibody Pair (STJAP00383)

STJAP00383

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Antibody Pairs  |
| <b>Short Description</b> | This ANGPTL3 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA. |
| <b>Applications</b>      | ELISA   |
| <b>Host/Source</b>       | Rat/Sheep   |
| <b>Reactivity</b>        | Human   |

### PRODUCT PROPERTIES

|                       |   |
|-----------------------|---|
| <b>Clonality</b>      |   |
| <b>Clone ID</b>       |   |
| <b>Concentration</b>  | 0.5 mg/mL   |
| <b>Conjugation</b>    | Unconjugated/Biotin   |
| <b>Purification</b>   | Capture Purification: Protein A or G.<br>Detection Purification: Antigen Affinity.  |
| <b>Dilution Range</b> |   |
| <b>Formulation</b>    | Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5.<br>Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5. |
| <b>Isotype</b>        |   |
| <b>Storage</b>        | Store at -20°C for long term, and avoid freeze-thaw cycles.   |
| <b>Instruction</b>    |   |

### TARGET INFORMATION

|                    |  |
|--------------------|--|
| <b>Gene ID</b>     | <a href="#">27329</a>  |
| <b>Gene Symbol</b> | <a href="#">ANGPTL3</a>  |
| <b>Uniprot ID</b>  | <a href="#">ANGL3_HUMAN</a>  |
| <b>Immunogen</b>   | Capture immunogen: Recombinant Human ANGPTL3 protien.<br>Detection Immunogen: Recombinant Human ANGPTL3 protien. |
| <b>Immunogen</b>   |  |
| <b>Region</b>      |  |
| <b>Specificity</b> |  |
| <b>Immunogen</b>   |  |
| <b>Sequence</b>    |  |