

Human IGFBP-2 (Insulin-Like Growth Factor-Binding Protein 2) Antibody Pair (STJAP00400)

STJAP00400

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This IGFBP-2 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3485
Gene Symbol	IGFBP2
Uniprot ID	IBP2_HUMAN
Immunogen	Capture immunogen: Recombinant Human IGFBP-2 protien. Detection Immunogen: Recombinant Human IGFBP-2 protien.
Immunogen Region	
Specificity	
Immunogen Sequence	