

Human THP (Uromodulin) Antibody Pair (STJAP00423)

STJAP00423

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This THP antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Sheep
Reactivity	Human

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7369
Gene Symbol	UMOD
Uniprot ID	UROM_HUMAN
Immunogen	Capture immunogen: Recombinant Human THP protien. Detection Immunogen: Recombinant Human THP protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	