

Human PD-1 (Programmed Cell Death Protein 1) Antibody Pair (STJAP00431)

STJAP00431

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This PD-1 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5133
Gene Symbol	PDCD1
Uniprot ID	PDCD1_HUMAN
Immunogen	Capture immunogen: Recombinant Human PD-1 protien. Detection Immunogen: Recombinant Human PD-1 protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	