

Human VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) Antibody Pair (STJAP00442)

STJAP00442

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This VEGFR-2/KDR antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rabbit/Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3791
Gene Symbol	KDR
Uniprot ID	VGFR2_HUMAN
Immunogen	Capture immunogen: Recombinant Human VEGFR-2/KDR protien. Detection Immunogen: Recombinant Human VEGFR-2/KDR protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	