

Human IL-1ra/IL-1F3 (Interleukin-1 Receptor Antagonist Protein) Antibody Pair (STJAP00448)

STJAP00448

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This IL-1ra/IL-1F3 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3557
Gene Symbol	IL1RN
Uniprot ID	IL1RA_HUMAN
Immunogen	Capture immunogen: Recombinant Human IL-1ra/IL-1F3 protien. Detection Immunogen: Recombinant Human IL-1ra/IL-1F3 protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	