

Human PPY (Pancreatic Polypeptide Prohormone) Antibody Pair (STJAP00451)

STJAP00451

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This PPY antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rabbit/Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5539
Gene Symbol	PPY
Uniprot ID	PAHO_HUMAN
Immunogen	Capture immunogen: Recombinant Human PPY protien. Detection Immunogen: Recombinant Human PPY protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	