

Human NEFL (Neurofilament Light Polypeptide) Antibody Pair (STJAP00455)

STJAP00455

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This NEFL antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rat/Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4747
Gene Symbol	NEFL
Uniprot ID	NFL_HUMAN
Immunogen	Capture immunogen: Natural Human NEFL protien. Detection Immunogen: Natural Human NEFL protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	