

Human C3a (Complement C3) Antibody Pair (STJAP00459)

STJAP00459

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This C3a antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	718
Gene Symbol	C3
Uniprot ID	CO3_HUMAN
Immunogen	Capture immunogen: Recombinant Human C3a protien. Detection Immunogen: Recombinant Human C3a protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	