

Human ISG15 (Ubiquitin-Like Protein Isg15) Antibody Pair (STJAP00483)

STJAP00483

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This ISG15 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9636
Gene Symbol	ISG15
Uniprot ID	ISG15_HUMAN
Immunogen	Capture immunogen: Recombinant Human ISG15 protien. Detection Immunogen: Recombinant Human ISG15 protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	