

Human PL (Chorionic Somatomammotropin Hormone 1) Antibody Pair (STJAP00488)

STJAP00488

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This PL antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Rabbit/Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1442
Gene Symbol	CSH1
Uniprot ID	CSH1_HUMAN
Immunogen	Capture immunogen: Recombinant Human PL protien. Detection Immunogen: Recombinant Human PL protien.
Immunogen Region	
Specificity	
Immunogen Sequence	