

Mouse IP-10/CXCL10 (C-X-C Motif Chemokine 10) Antibody Pair (STJAP00522)

STJAP00522

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This IP-10/CXCL10 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	15945
Gene Symbol	Cxcl10
Uniprot ID	CXL10_MOUSE
Immunogen	Capture immunogen: Recombinant Mouse IP-10/CXCL10 protien. Detection Immunogen: Recombinant Mouse IP-10/CXCL10 protien.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	